

VERY IMPORTANT

Cleaning Procedures

Normal cleaning procedures:

May be cleaned with mild dish soap and warm to hot water. Use only a soft cloth to prevent scratching. Do not use abrasive cleaners or scouring pads.

Disinfecting cleaning procedures:

Use 1:10 chlorine bleach to water solution (13 ounces bleach in 1 gallon water). Many manufacturers produce various types of products with different chemicals. Not all chemicals have been tested against our product. Testing is ongoing and if a chemical is found to have an adverse reaction with our product, we will place it on the Do Not Use list.

There is no specific product that is recommended or endorsed by Paramedic Design, Inc.

Do Not Use any product containing the chemical:

Ammonia or Ammonium Chloride - known to cause cracking/crazing

PARAMEDIC DESIGN, INC.
— Acrylic EMS Organizers — 

Phone: 1-800-527-5943

Web: www.paramedicdesigninc.com

Email: products@paramedicdesigninc.com